

Hill-Rom.

Affinity® Birthing Bed
Tempur-Pedic® Mattress



You've created the most comfortable environment.
We've created the most comfortable mattress.

Tempur-Pedic's sleep-science technology—ideal comfort to improve expectant and postpartum moms' sleep and health.

Our TEMPUR® material products deliver a relaxing, rejuvenating quality of sleep every new mom needs!

Better sleep and greater comfort equate to greater patient satisfaction.

- 64% of patients said the bed makes them feel like the hospital cares about them.
- 81% of patients said the mattress helps them to get comfortable and relax in the hospital environment.
- 83% of patients said the mattress helped them to sleep soundly in spite of noise and interruptions in the hospital environment.

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom®

Affinity® Birthing Bed Tempur-Pedic® Mattress



Technical Specifications

Affinity® Birthing Bed Tempur-Pedic® Mattress

Multi-Layer Foam Construction

Dual-Zone Top Layer:

Upper Body Soft Support; Density 5.3 PCF

Upper Body Base 29 IFD, Density 1.85 PCF

Seat Area:

Transitions into a triple laminated core. This section consists of a transitional wedge (3" resolving to 1" at the end of the mattress) of T85-13. This material is positioned over a 2" layer of TEMPUR T110-245, a 7 PCF material used in the manufacture of high performance wheelchair cushions. The reverse wedge (1" - 3") in the mattress substrate is manufactured using Ensolite foam.

Upper Body:

Length 51.5"

Width 35.25" +/- .5"

Height 5"

Weight Approx. 20 lbs.

Weight Limit 500 lbs.

Flammability Cal 117 standard, call 129, 16 CFR 1633

Warranty 3 years

MVP 10 gm/ml2/hr

Foot Section is not designed with Tempur-Pedic material.

Non-Powered Model

3" layer of TEMPUR® material in Upper Torso Area.

Triple laminated core in seat area including 2" of high density TEMPUR® material used in premium wheelchair cushions.

Dual cover design. Both the inner and outer covers are of a specially developed material from Uretek. Seams sewn and RF welded.

Benefit

The material used in Tempur-Pedic® medical mattresses is **clinically proven to reduce the incidence of stage II or higher ulcers. Surveys show patients are more comfortable and feel the hospital really cares about their satisfaction.**

Provides **support, stability**, and maintains **superior levels of pressure redistribution.**

- **Superior hydrolysis resistance**
- **Withstands harsh cleaning exceptionally well**
- **Seams flexible and fully sealed**

Air Assist Model

3" layer of TEMPUR® material in Upper Torso Area.

Triple laminated core in seat area including 2" of high density TEMPUR® material used in premium wheelchair cushions.

Unique one-piece air bladder.

Dual cover design. Both the inner and outer covers are of a specially developed material from Uretek. Seams sewn and RF welded.

Benefit

The material used in Tempur-Pedic® medical mattresses is **clinically proven to reduce the incidence of stage II or higher ulcers. Surveys show patients are more comfortable and feel the hospital really cares about their satisfaction.**

Provides **support, stability**, and maintains **superior levels of pressure redistribution.**

Provides **lumbar, elevation and seated support** as well as **support for positioning.**

- **Superior hydrolysis resistance**
- **Withstands harsh cleaning exceptionally well**
- **Seams flexible and fully sealed**

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

152814 rev 2 5/22/08

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

USA **800-445-3730**
Canada **800-267-2337**
www.hill-rom.com